



Firefly Luciferase + eGFP Lentifact™ Purified Lentiviral Particles • Cat No. LPP-HLUC-LV201-025

Ready-to-use lentiviral particles for the transduction of a variety of mammalian cells including difficult-to-transfect, primary, stem and non-dividing cells.

Description

- Produced with a standardized protocol using highly purified plasmids and EndoFectin Lenti™ transfection and TiterBoost™ reagents
- CMV promoter for the expression of hLUC (humanized firefly luciferase) in target cells
- SV40 promoter for the expression of GFP
- Puromycin resistance marker (*Pac* gene) for selection of stably transduced cells
- *Pac* gene bicistronically coexpressed with GFP

Contents and storage

Provided as 1 vial of 25 µl of purified lentiviral particles with titers of $\sim 3 \times 10^9$ copies/ml.

Lentifact particles are shipped on dry ice and **must be stored at -80°C immediately upon receipt**. Avoid repeated freeze-thaw cycles as this will reduce titers.

Quality control

The lentiviral expression construct was validated by full-length sequencing, restriction enzyme digestion and PCR-size validation using gene-specific and vector-specific primers. Product is confirmed free of bacteria, fungi and common *Mycoplasma* contamination.

Viral titer

Lentivirus products were titrated using qRT-PCR, which determines the physical copy numbers of viral genomic RNA. We suggest that the customer estimate the transduction unit (TU or IFU) for the specific host cells based on the following formula before transduction:

$\text{TU} = \text{Titer (physical copy number)} / 100$.

The customer should test the transduction at MOI=0.3, 1, 3, 5, 10 using TU for the best transduction efficiency.

Overview of production

The OmicsLink™ ORF lentiviral expression plasmid (GeneCopoeia Cat. No. EX-HLUC-LV201) was constructed using GeneCopoeia proprietary RecJoin™ technology. This plasmid was co-transfected into 293Ta cells (GeneCopoeia Cat. No. CLV-PK-01) with the Lenti-Pac HIV Packaging Mix (GeneCopoeia Cat. No. HPK-LvTR-20). Lentivirus-containing supernatants were harvested 48 hours after transfection, the particles were purified and stored at -80°C in aliquots.

User manual

Please contact GeneCopoeia for a copy or download at:
http://genecopoeia.com/product/lentiviral/pdf/packaging_kit_manual.pdf

GeneCopoeia, Inc.
9620 Medical Center Drive, #101
Rockville, Maryland 20850
Tel: 301-762-0888 Fax: 301-762-8333
Email: inquiry@genecopoeia.com
Web: www.genecopoeia.com

GeneCopoeia Products are for Research Use Only
Trademarks: Lentifact™, Lenti-Pac™, OmicsLink™, EndoFectin™, TiterBoost™ (GeneCopoeia Inc.)

Copyright © 2011 GeneCopoeia Inc.
LPPHLUCGFP090211