

Plate	Position	Catalog No. of Primer	Accession No. of Gene	Symbol
QG056-01	A01	HQP009024	NM_000038	APC
QG056-01	A02	HQP006794	NM_015399	BRMS1
QG056-01	A03	HQP016627	NM_006273	CCL7
QG056-01	A04	HQP022972	NM_000610	CD44
QG056-01	A05	HQP023466	NM_004360	CDH1
QG056-01	A06	HQP000131	NM_001797	CDH11
QG056-01	A07	HQP000065	NM_004932	CDH6
QG056-01	A08	HQP000369	NM_000077	CDKN2A
QG056-01	A09	HQP001249	NM_001273	CHD4
QG056-01	A10	HQP002506	NM_001846	COL4A2
QG056-01	A11	HQP021171	NM_003650	CST7
QG056-01	A12	HQP003515	NM_001903	CTNNA1
QG056-01	B01	HQP003626	NM_000396	CTSK
QG056-01	B02	HQP003630	NM_001912	CTSL1
QG056-01	B03	HQP016669	NM_000609	CXCL12
QG056-01	B04	HQP018803	NM_003467	CXCR4
QG056-01	B05	HQP021285	NM_003677	DENR
QG056-01	B06	HQP004939	NM_004442	EPHB2
QG056-01	B07	HQP053982	NM_001986	ETV4
QG056-01	B08	HQP005034	NM_005243	EWSR1
QG056-01	B09	HQP005141	NM_005245	FAT
QG056-01	B10	HQP005439	NM_002011	FGFR4
QG056-01	B11	HQP005915	NM_002020	FLT4
QG056-01	B12	HQP006022	NM_002026	FN1
QG056-01	C01	HQP013317	NM_014164	FXYD5
QG056-01	C02	HQP007777	NM_000825	GNRH1
QG056-01	C03	HQP020835	NM_032551	KISS1R
QG056-01	C04	HQP008800	NM_000601	HGF
QG056-01	C05	HQP000974	NM_006665	HPSE
QG056-01	C06	HQP009036	NM_005343	HRAS
QG056-01	C07	HQP000660	NM_006410	HTATIP2
QG056-01	C08	HQP009518	NM_000618	IGF1
QG056-01	C09	HQP009718	NM_001562	IL18
QG056-01	C10	HQP009641	NM_000576	IL1B
QG056-01	C11	HQP009681	NM_001557	IL8RB
QG056-01	C12	HQP009800	NM_002206	ITGA7
QG056-01	D01	HQP009818	NM_000212	ITGB3
QG056-01	D02	HQP009861	NM_002231	CD82
QG056-01	D03	HQP010098	NM_002256	KISS1
QG056-01	D04	HQP053983	NM_002295	RPSA
QG056-01	D05	HQP011099	NM_006500	MCAM

Plate	Position	Catalog No. of Primer	Accession No. of Gene	Symbol
QG056-01	D06	HQP011135	NM_002392	MDM2
QG056-01	D07	HQP011181	NM_000245	MET
QG056-01	D08	HQP001123	NM_006838	METAP2
QG056-01	D09	HQP011196	NM_002410	MGAT5
QG056-01	D10	HQP011264	NM_002425	MMP10
QG056-01	D11	HQP011265	NM_005940	MMP11
QG056-01	D12	HQP011257	NM_002422	MMP3
QG056-01	E01	HQP011258	NM_002423	MMP7
QG056-01	E02	HQP011263	NM_004994	MMP9
QG056-01	E03	HQP022096	NM_004689	MTA1
QG056-01	E04	HQP023201	NM_014751	MTSS1
QG056-01	E05	HQP011597	NM_002467	MYC
QG056-01	E06	HQP011601	NM_005376	MYCL1
QG056-01	E07	HQP011778	NM_000268	NF2
QG056-01	E08	HQP011853	NM_000269	NME1
QG056-01	E09	HQP053986	NM_002512	NME2
QG056-01	E10	HQP011857	NM_005009	NME4
QG056-01	E11	HQP019536	NM_006981	NR4A3
QG056-01	E12	HQP013207	NM_002659	PLAUR
QG056-01	F01	HQP013392	NM_002687	PNN
QG056-01	F02	HQP015535	NM_000314	PTEN
QG056-01	F03	HQP016131	NM_000321	RB1
QG056-01	F04	HQP016377	NM_006914	RORB
QG056-01	F05	HQP016841	NM_003011	SET
QG056-01	F06	HQP054007	NM_005901	SMAD2
QG056-01	F07	HQP010961	NM_005359	SMAD4
QG056-01	F08	HQP017696	NM_005417	SRC
QG056-01	F09	HQP017744	NM_001050	SSTR2
QG056-01	F10	HQP017845	NM_003177	SYK
QG056-01	F11	HQP017971	NM_005650	TCF20
QG056-01	F12	HQP018044	NM_000660	TGFB1
QG056-01	G01	HQP018093	NM_003255	TIMP2
QG056-01	G02	HQP018094	NM_000362	TIMP3
QG056-01	G03	HQP018095	NM_003256	TIMP4
QG056-01	G04	HQP021502	NM_003810	TNFSF10
QG056-01	G05	HQP018175	NM_000546	TP53
QG056-01	G06	HQP011252	NM_002420	TRPM1
QG056-01	G07	HQP018271	NM_000369	TSHR
QG056-01	G08	HQP018481	NM_003376	VEGFA
QG056-01	G09	HQP022974	NM_001001390	CD44
QG056-01	G10	HQP009026	NM_000194	HPRT1

Plate	Position	Catalog No. of Primer	Accession No. of Gene	Symbol
QG056-01	G11	HQP006023	NM_054034	FN1
QG056-01	G12	HQP018272	NM_001018036	TSHR
QG056-01	H01	HGDC		
QG056-01	H02	HGDC		
QG056-01	H03	HQP006940	NM_002046	GAPDH
QG056-01	H04	HQP016381	NM_001101	ACTB
QG056-01	H05	HQP015171	NM_004048	B2M
QG056-01	H06	HQP006171	NM_012423	RPL13A
QG056-01	H07	HQP009026	NM_000194	HPRT1
QG056-01	H08	HQP054253	NR_003286	RN18S1
QG056-01	H09	RT		
QG056-01	H10	RT		
QG056-01	H11	PCR		
QG056-01	H12	PCR		