

Plate	Position	Catalog No.	Mature miRNA Accession No.	Mature miRNA ID	miRNAs also Amplified
QM018-01	A01	HmiRQP0044	MIMAT0000416	hsa-miR-1	
QM018-01	A02	HmiRQP0021	MIMAT0000099	hsa-miR-101-3p	
QM018-01	A03	HmiRQP0096	MIMAT0000423	hsa-miR-125b-5p	
QM018-01	A04	HmiRQP0192	MIMAT0000437	hsa-miR-145-5p	
QM018-01	A05	HmiRQP0354	MIMAT0000082	hsa-miR-26a-5p	
QM018-01	A06	HmiRQP0396	MIMAT0000244	hsa-miR-30c-5p	
QM018-01	A07	HmiRQP0137	MIMAT0000242	hsa-miR-129-5p	
QM018-01	A08	HmiRQP0316	MIMAT0000076	hsa-miR-21-5p	
QM018-01	A09	HmiRQP0317	MIMAT0000267	hsa-miR-210	
QM018-01	A10	HmiRQP0349	MIMAT0000080	hsa-miR-24-3p	
QM018-01	A11	HmiRQP0099	MIMAT0000445	hsa-miR-126-3p	
QM018-01	A12	HmiRQP0188	MIMAT0000435	hsa-miR-143-3p	
QM018-01	B01	HmiRQP0247	MIMAT0000455	hsa-miR-185-5p	
QM018-01	B02	HmiRQP0283	MIMAT0000461	hsa-miR-195-5p	
QM018-01	B03	HmiRQP0307	MIMAT0000266	hsa-miR-205-5p	
QM018-01	B04	HmiRQP0338	MIMAT0000278	hsa-miR-221-3p	
QM018-01	B05	HmiRQP0342	MIMAT0000280	hsa-miR-223-3p	
QM018-01	B06	HmiRQP0002	MIMAT0000062	hsa-let-7a-5p	
QM018-01	B07	HmiRQP0012	MIMAT0000067	hsa-let-7f-5p	
QM018-01	B08	HmiRQP0510	MIMAT0001635	hsa-miR-452-5p	
QM018-01	B09	HmiRQP0018	MIMAT0000098	hsa-miR-100-5p	
QM018-01	B10	HmiRQP0509	MIMAT0001631	hsa-miR-451a	
QM018-01	B11	HmiRQP0854	MIMAT0000097	hsa-miR-99a-5p	
QM018-01	B12	HmiRQP0400	MIMAT0000692	hsa-miR-30e-5p	
QM018-01	C01	HmiRQP0401	MIMAT0000089	hsa-miR-31-5p	
QM018-01	C02	HmiRQP0357	MIMAT0000083	hsa-miR-26b-5p	
QM018-01	C03	HmiRQP0497	MIMAT0001536	hsa-miR-429	
QM018-01	C04	HmiRQP0786	MIMAT0000252	hsa-miR-7-5p	
QM018-01	C05	HmiRQP0305	MIMAT0000264	hsa-miR-203	
QM018-01	C06	HmiRQP0339	MIMAT0000279	hsa-miR-222-3p	
QM018-01	C07	HmiRQP0221	MIMAT0000646	hsa-miR-155-5p	
QM018-01	C08	HmiRQP0110	MIMAT0004604	hsa-miR-127-5p	
QM018-01	C09	HmiRQP0346	MIMAT0000418	hsa-miR-23b-3p	
QM018-01	C10	HmiRQP0139	MIMAT0004605	hsa-miR-129-2-3p	
QM018-01	C11	HmiRQP0111	MIMAT0000446	hsa-miR-127-3p	
QM018-01	C12	HmiRQP0007	MIMAT0000065	hsa-let-7d-5p	
QM018-01	D01	HmiRQP0298	MIMAT0000682	hsa-miR-200a-3p	
QM018-01	D02	HmiRQP0184	MIMAT0000432	hsa-miR-141-3p	
QM018-01	D03	HmiRQP0300	MIMAT0000318	hsa-miR-200b-3p	
QM018-01	D04	HmiRQP0302	MIMAT0000617	hsa-miR-200c-3p	
QM018-01	D05	HmiRQP0032	MIMAT0000253	hsa-miR-10a-5p	

Plate	Position	Catalog No.	Mature miRNA Accession No.	Mature miRNA ID	miRNAs also Amplified
QM018-01	D06	HmiRQP0575	MIMAT0002852	hsa-miR-517a-3p	
QM018-01	D07	HmiRQP0230	MIMAT0000070	hsa-miR-17-5p	
QM018-01	D08	HmiRQP0375	MIMAT0000681	hsa-miR-29c-3p	
QM018-01	D09	HmiRQP0361	MIMAT0000419	hsa-miR-27b-3p	
QM018-01	D10	HmiRQP0022	MIMAT0000101	hsa-miR-103a-3p	
QM018-01	D11	HmiRQP0167	MIMAT0000770	hsa-miR-133b	
QM018-01	D12	HmiRQP0344	MIMAT0000078	hsa-miR-23a-3p	
QM018-01	E01	HmiRQP0166	MIMAT0000427	hsa-miR-133a	
QM018-01	E02	HmiRQP0289	MIMAT0000232	hsa-miR-199a-3p	
QM018-01	E03	HmiRQP0337	MIMAT0004568	hsa-miR-221-5p	
QM018-01	E04	HmiRQP0299	MIMAT0004571	hsa-miR-200b-5p	
QM018-01	E05	HmiRQP0001	MIMAT0004481	hsa-let-7a-3p	
QM018-01	E06	HmiRQP0008	MIMAT0004484	hsa-let-7d-3p	
QM018-01	E07	HmiRQP0011	MIMAT0004486	hsa-let-7f-1-3p	
QM018-01	E08	HmiRQP0019	MIMAT0004512	hsa-miR-100-3p	
QM018-01	E09	HmiRQP0020	MIMAT0004513	hsa-miR-101-5p	
QM018-01	E10	HmiRQP0031	MIMAT0004555	hsa-miR-10a-3p	
QM018-01	E11	HmiRQP0136	MIMAT0004548	hsa-miR-129-1-3p	
QM018-01	E12	HmiRQP0187	MIMAT0004599	hsa-miR-143-5p	
QM018-01	F01	HmiRQP0191	MIMAT0004601	hsa-miR-145-3p	
QM018-01	F02	HmiRQP0220	MIMAT0004658	hsa-miR-155-3p	
QM018-01	F03	HmiRQP0301	MIMAT0004657	hsa-miR-200c-5p	
QM018-01	F04	HmiRQP0341	MIMAT0004570	hsa-miR-223-5p	
QM018-01	F05	HmiRQP0345	MIMAT0004496	hsa-miR-23a-5p	
QM018-01	F06	HmiRQP0347	MIMAT0004587	hsa-miR-23b-5p	
QM018-01	F07	HmiRQP0348	MIMAT0000079	hsa-miR-24-1-5p	
QM018-01	F08	HmiRQP0360	MIMAT0004588	hsa-miR-27b-5p	
QM018-01	F09	HmiRQP0376	MIMAT0004673	hsa-miR-29c-5p	
QM018-01	F10	HmiRQP0402	MIMAT0004504	hsa-miR-31-3p	
QM018-01	F11	HmiRQP0576	MIMAT0002851	hsa-miR-517-5p	
QM018-01	F12	HmiRQP0855	MIMAT0004511	hsa-miR-99a-3p	
QM018-01	G01	HmiRQP0095	MIMAT0004592	hsa-miR-125b-1-3p	
QM018-01	G02	HmiRQP0183	MIMAT0004598	hsa-miR-141-5p	
QM018-01	G03	HmiRQP0315	MIMAT0004494	hsa-miR-21-3p	
QM018-01	G04	HmiRQP0340	MIMAT0004569	hsa-miR-222-5p	
QM018-01	G05	HmiRQP0356	MIMAT0004500	hsa-miR-26b-3p	
QM018-01	G06	HmiRQP0390	MIMAT0000088	hsa-miR-30a-3p	
QM018-01	G07	HmiRQP0511	MIMAT0001636	hsa-miR-452-3p	
QM018-01	G08				
QM018-01	G09				
QM018-01	G10				

Plate	Position	Catalog No.	Mature miRNA Accession No.	Mature miRNA ID	miRNAs also Amplified
QM018-01	G11				
QM018-01	G12				
QM018-01	H01	NC	NC	NC	
QM018-01	H02	NC	NC	NC	
QM018-01	H03	HmiRQP9001	NR_002752	HK1	
QM018-01	H04	HmiRQP9011	NR_002750	HK2	
QM018-01	H05	HmiRQP9021	NR_002745	HK3	
QM018-01	H06	HmiRQP9051	NR_002746	HK4	
QM018-01	H07	HmiRQP9061	NR_002744	HK5	
QM018-01	H08	HmiRQP9071	NR_002450	HK6	
QM018-01	H09	RT	RT	RT	
QM018-01	H10	RT	RT	RT	
QM018-01	H11	PCR	PCR	PCR	
QM018-01	H12	PCR	PCR	PCR	